

**NEW**




# NECK CONTOUR PACK

ENCIRCLES  
NECK  
COMPLETELY



COVERS  
ENTIRE  
CERVICAL  
AREA

CAN BE  
APPLIED  
DURING  
TRACTION



**A HYDROCOLLATOR® Steam Pack specially designed for convenient, effective application to the neck.**

**A tremendous aid in the treatment of neck injuries and general reduction of spasm and pain in the neck area.**

**Shaped to cover the entire cervical area, and to conform comfortably around the entire neck.**

**Used with the special Hydrocollator Cover-Sling, this Pack is conveniently held in place -- regardless of the patient's position. Particularly helpful with the patient in the supine position -- usual position in neck treatment.**

**Also, with the special cover, the Neck Contour Pack can be applied during traction -- a valuable adjunct to cervical traction. Users report definitely enhanced comfort of their patients during traction and improved quality of their cervical traction results.**

**Hydrocollator quality throughout.**

U. S. PAT. No. 2,710,008

D-424

## HYDROCOLLATOR STEAM PACKS

SIZE	DIMENSIONS (DRY)	SHAPE	PRICE FOB CHATTANOOGA
STANDARD SIZE	10" x 12"		\$4.95 Ea.
HALF SIZE	5" x 12"		\$2.75 Ea.
OVERSIZE	24" x 15"		\$10.50 Ea.
NECK CONTOUR	24" Long		\$6.00 Ea.
EYE SIZE	3" x 3 1/2"		\$1.85 BOX OF TWO

**OVER**  
FOR DETAILS  
ON THE  
NEW  
COVER  
SLING

ORIGINATED AND MFD. BY  
**CHATTANOOGA PHARMACAL CO., INC.**  
CHATTANOOGA 15, TENNESSEE




®

HYDROCOLLATOR


# COVER-SLING

FOR NECK CONTOUR PACK



**Apply your Neck Contour Packs more conveniently and effectively with this efficient Terry Cloth Cover-Sling.**

**With far less unwieldy bulk and mass than the usual multiple towel layers, the special foam filler gives much closer conformation to the neck area. This foam filler provides the proper flesh-to-heated-pack spacing, and its porosity is just right to permit penetration of the steam from the heated Pack (for true Moist Heat applications).**



**The double pocket design holds the foam filler in one pocket and the heated Neck Contour Pack in the other. When folded over, you have a neat, convenient-to-apply, wrap-around Pack.**

**Easily adjustable to any neck size.**

**Regardless of the patient's position or movements, the Pack is held firmly in place, thus, effective Hydrocollator® treatment can be continued during cervical traction. The patient is free to move around with complete ease - while wearing the heated Neck Contour Pack.**

**Hydrocollator quality throughout.**

**PRICES:**

F. O. B. CHATTANOOGA

EACH

**\$2.50**

ORIGINATED AND MFD. BY  
**CHATTANOOGA PHARMACAL CO., INC.**  
 CHATTANOOGA 15, TENNESSEE